TESTS

REF 06454011

[Item name] Blood glucose test for personal use
[Product name] Accu-Chek Performa Strips
[Approval Number] in-vitro수허13-1083호

Suitable for self-testing

Intended Use

The Accu-Chek Performa test strips are intended to be used with the Accu-Chek Performa, Accu-Chek Performa Nano, Accu-Chek Performa Connect, and Accu-Chek Performa Solo blood glucose meters to quantitatively measure glucose in fresh capillary whole blood as an aid in monitoring the effectiveness of controlling glucose.

Testing sites for the Accu-Chek Performa, Accu-Chek Performa Nano, and Accu-Chek Performa Connect meters include the finger, palm, forearm, and upper arm. The testing site for the Accu-Chek Performa Solo meter is the fingertip.

The Accu-Chek Performa test strips, used with the Accu-Chek Performa, Accu-Chek Performa Nano, Accu-Chek Performa Connect, and Accu-Chek Performa Solo blood glucose meters, are intended for in vitro diagnostic use by people with diabetes at home. The Accu-Chek Performa Solo meter is for home use only.

The systems are not for use in diagnosis or screening of diabetes mellitus, nor for testing neonate cord blood samples.

Consumer Information

Important information: These test strips are labeled with a green symbol to distinguish them from earlier test strips that were subject to a clinically relevant maltose interference.* The green symbol can be found on the test strip box and on the label of the test strip container.

*Data on file

Warning: Choking hazard. Small parts. Keep away from children under the age of 3 years.

Warning: Inspect the test strip container before using the test strips for the first time. If you see any damage to the container, if anything prevents the cap from closing properly, or if the container was open before using for the first time, do not use the test strips. Do not perform a control test. Contact Roche. Damaged test strips can cause inaccurate results, which could lead to improper therapy.

Introduction

Testing your blood glucose regularly may help you better manage your diabetes. Medical studies show that you and your healthcare professional can manage your blood glucose to near normal levels.¹ This can prevent or slow the development of complications from diabetes.

Contents of the pack

Pack containing test strips and package inserts.

All components of the pack can be discarded in domestic waste. Because the reactive substances are in such small quantities, they are not considered to be hazardous materials under EU regulations. Discard used test strips according to local regulations. If you have any questions, contact Roche.

Test strip storage and handling

- Store the test strips at temperatures between 2 and 30 °C. Do not freeze the test strips.
- Use the test strips at temperatures between 8 and 44 °C.
- Use the test strips between 10 and 90 % humidity. Do not store the test strips in high heat and moisture areas such as the bathroom or kitchen.
- Store the unused test strips in their original test strip container with the cap closed.
- Close the test strip container tightly immediately after removing a test strip to protect the test strips from humidity.
- Use the test strip immediately after removing it from the test strip container.
- Discard the test strips if they are past the use by date. Expired test strips can produce incorrect results. The use by date is

printed on the test strip box and on the label of the test strip container next to . The test strips can be used until the printed use by date when they are stored and used correctly. This applies for test strips from a new, unopened test strip container and for test strips from a test strip container that has already been opened.

Performing a Blood Glucose Test

Note: If your meter requires an activation chip, contact Roche to obtain one.

If you have poor circulation, testing your own blood glucose may not be right for you. Ask your healthcare professional.

- Wash your hands in warm, soapy water. Rinse and dry completely.
- Prepare the lancing device.
- Check the use by date on the test strip container. Do not use test strips past the use by date.
- Insert the test strip into the meter in the direction of the arrows. The meter turns on.
- Obtain a blood drop using the lancing device.
- Touch the blood drop to the **front edge** of the yellow window of the test strip. Do not put blood on top of the test strip. When you see flash, there is enough blood in the test strip.

Understanding Test Results

The normal fasting glucose level for a non-diabetic adult is below 100 mg/dL (5.6 mmol/L).² The normal glucose level for a non-diabetic adult 2 hours post-meal, e.g. simulated by 75 g Oral Glucose Tolerance Test (OGTT), is less than 140 mg/dL (7.8 mmol/L).³ A criterion for the diagnosis of diabetes in adults is a fasting glucose level of 126 mg/dL or higher (7.0 mmol/L or higher) confirmed in two tests.^{2,4,5} Adults with a fasting glucose level between 100 and 125 mg/dL (5.6 and 6.9 mmol/L) are defined as having impaired fasting glucose (prediabetes).² Other diagnostic criteria for diabetes exist. Consult your healthcare professional to determine if you have diabetes or not. For people with diabetes: Consult your healthcare professional for the blood glucose range appropriate for you. You should treat your low or high blood glucose as recommended by your healthcare professional.

These test strips deliver results that correspond to blood glucose concentrations in plasma as per the recommendation of the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC).⁶ Therefore, the meter displays blood glucose concentrations that refer to plasma although whole blood is always applied to the test strip.

Unusual test results

If **LO** is displayed on the meter, your blood glucose may be below 10 mg/dL (0.6 mmol/L).

If **HI** is displayed on the meter, your blood glucose may be over 600 mg/dL (33.3 mmol/L).

For detailed information on error messages, refer to the User's Manual.

If your blood glucose result does not match how you feel, follow these steps:

- Repeat the blood glucose test with a new test strip.
- Perform a control test with an Accu-Chek Performa control solution.
- Check this list to help solve the problem.
 - Were the test strips expired?
 - Was the cap on the test strip container always closed tightly?
 - Was the test strip used immediately after it was removed from the test strip container?
 - Were the test strips stored in a cool, dry place?
 - Did you follow the directions?
- If you think your blood glucose results are too low, too high, or doubtful, contact your healthcare professional.

Limitations

- Blood concentrations of galactose >15 mg/dL (>0.83 mmol/L) will cause overestimation of blood glucose results.
- Lipemic samples (triglycerides) >1,800 mg/dL (>20.3 mmol/L) may produce elevated blood glucose results.
- Intravenous administration of ascorbic acid which results in blood concentrations of ascorbic acid >3 mg/dL (>0.17 mmol/L) will cause overestimation of blood glucose results.

- Intravenous administration of N-acetylcysteine, which results in blood concentrations of N-acetylcysteine >5 mg/dL (>307 µmol/L), will cause overestimation of blood glucose results. Do not use during intravenous administration of high-dose N-acetylcysteine.
- If peripheral circulation is impaired, collection of capillary blood from the approved sample sites is not advised as the results might not be a true reflection of the physiological blood glucose level. This may apply in the following circumstances: Severe dehydration as a result of diabetic ketoacidosis or due to hyperglycemic hyperosmolar non-ketotic syndrome, hypotension, shock, decompensated heart failure NYHA Class IV, or peripheral arterial occlusive disease.
- Hematocrit should be between 10 and 65 %. Ask your healthcare professional if you do not know your hematocrit.
- This system has been tested at altitudes up to 3,094 meters.

Performance Characteristics

The Accu-Chek Performa system complies with the requirements of ISO 15197:2013 (In vitro diagnostic test systems – Requirements for blood glucose monitoring systems for self-testing in managing diabetes mellitus).

Calibration and traceability: The system (meter and test strips) is calibrated with venous blood containing various glucose concentrations as a calibrator. The reference values are obtained using the hexokinase method which is calibrated using the ID-GCMS method. The ID-GCMS method as the method of highest metrological quality (order) is traceable to a primary NIST standard. Using this traceability chain, the results obtained with these test strips for control solutions can also be traced back to the NIST standard.

Detection limit (lowest value displayed): 10 mg/dL (0.6 mmol/L) for the test strip

System measurement range: 10–600 mg/dL (0.6–33.3 mmol/L)

Sample size: 0.6 µL

Test time: 5 seconds

System accuracy:

System accuracy results for glucose concentrations less than 100 mg/dL (less than 5.55 mmol/L)

within ±5 mg/dL (within ±0.28 mmol/L)	within ±10 mg/dL (within ±0.56 mmol/L)	within ±15 mg/dL (within ±0.83 mmol/L)
137/168 (81.5 %)	163/168 (97.0 %)	167/168 (99.4 %)

System accuracy results for glucose concentrations equal to or greater than 100 mg/dL (equal to or greater than 5.55 mmol/L)

within ±5 %	within ±10 %	within ±15 %
256/432 (59.3 %)	393/432 (91.0 %)	428/432 (99.1 %)

System accuracy results for glucose concentrations between 12 mg/dL (0.67 mmol/L) and 525 mg/dL (29.1 mmol/L)

within ±15 mg/dL or within ±15 % (within ±0.83 mmol/L or within ±15 %)
595/600 (99.2 %)

Repeatability:

Mean value	[mg/dL]	42.6	89.5	122.4	186.9	313.0
	[mmol/L]	2.4	5.0	6.8	10.4	17.4
Standard deviation	[mg/dL]	1.7	3.2	4.3	7.2	11.5
	[mmol/L]	0.1	0.2	0.2	0.4	0.6
Coefficient of variation [%]		—	—	3.5	3.9	3.7

Intermediate precision:

Mean value	[mg/dL]	45.6	118.6	310.6
	[mmol/L]	2.5	6.6	17.2
Standard deviation	[mg/dL]	1.3	2.7	6.3
	[mmol/L]	0.09	0.2	0.3
Coefficient of variation [%]		—	2.3	2.0

Performance assessment by the user: A study evaluating glucose values from fingertip capillary blood samples obtained by 209 lay persons showed the following results:

- For glucose concentrations less than 100 mg/dL (less than 5.55 mmol/L), 97.6 % of the test results were within ±15 mg/dL (within ±0.83 mmol/L) of the results obtained through laboratory testing.

- For glucose concentrations equal to or greater than 100 mg/dL (equal to or greater than 5.55 mmol/L), 97.0 % of the test results were within ±15 % of the results obtained through laboratory testing.

Test principle: The enzyme on the test strip, mutant variant of quinoprotein glucose dehydrogenase (Mut. Q-GDH) from *Acinetobacter calcoaceticus*, recombinant in *E. coli*, converts the glucose in the blood sample to gluconolactone. This reaction creates a harmless DC electrical current that the meter interprets for the blood glucose result. The sample and environmental conditions are evaluated using AC and DC signals.

Reagent composition

Mediator	6.72 %
Quinoprotein glucose dehydrogenase	15.27 %
Pyrroloquinoline quinone	0.14 %
Buffer	34.66 %
Stabilizer	0.54 %
Non-reactive ingredients	42.66 %

Minimum at time of manufacture

From *A. calcoaceticus*, recombinant in *E. coli*, detailed description in patent application WO 2007/118647 (as “mutant 31” in table 4)

Note: For an explanation of symbols used and a list of references, refer to the end of this package insert.

Control and linearity test kits (if available)

Accu-Chek Performa control solution – Refer to the control solution package insert for details.

Accu-Chek linearity test kit – Refer to the linearity test kit package insert for details.

Visit our website at www.accu-chek.com or contact the local Roche representative for more information.

LAST UPDATE: 2022-02

References

- Davies MJ, D'Alessio DA, Fradkin J, Kernan WN, Mathieu C, Mingrone G, Rossing P, Tsapas A, Wexler DJ, Buse JB. Management of Hyperglycemia in Type 2 Diabetes, 2018. A Consensus Report by the American Diabetes Association (ADA) and the European Association for the Study of Diabetes (EASD). *Diabetes Care*. 2018 Dec; 41(12):2669-2701.
- American Diabetes Association: 2. Classification and Diagnosis of Diabetes: Standards of Medical Care in Diabetes-2019. *Diabetes Care*, 42, (Suppl. 1): S13-S28 (2019).
- American Diabetes Association website: Diagnosing Diabetes and Learning about Prediabetes. <http://www.diabetes.org/diabetes-basics/diagnosis/>. Accessed April 22, 2019.
- IDF Clinical Guidelines Task Force. Global guideline for Type 2 diabetes. Brussels: International Diabetes Federation, 2012.
- Definition and diagnosis of diabetes mellitus and intermediate hyperglycemia: report of a WHO/IDF consultation. WHO, Geneva 2006 (ISBN 92 4 159493 4, ISBN 978 92 4 159493 6).
- D'Orazio et al.: “Approved IFCC Recommendation on Reporting Results for Blood Glucose (Abbreviated);” *Clinical Chemistry* 51:9 1573-1576 (2005).

	Consult package insert
	Temperature limitation (store at)
	Use by (opened / unopened)
	Manufacturer
	Catalogue number
	Batch code
	In vitro diagnostic medical device
	Global Trade Item Number
	Serial number
	These test strips deliver results that refer to plasma as per IFCC, and the symbol distinguishes them from earlier test strips that were subject to a clinically relevant maltose interference.
	All components of the pack can be discarded in domestic waste. Discard used test strips according to local regulations.

ACCUCHEK Performa

08985111061A-0322

ACCUCHEK Performa

08985111061A-0322

Roche Diabetes Care GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany
www.accu-chek.com

IVD

ACCUCHEK Performa

08985111061A-0322

Roche

ACCUCHEK, ACCUCHEK PERFORMA, PERFORMA CONNECT, and PERFORMA NANO are trademarks of Roche.

© 2022 Roche Diabetes Care

This file may not print or view at 100%.

Color breaks do not print.

Roche USA – 84395
V1/1 – 08985111061A – Black

Far right column of front side and far left column of back side are keep out areas needed for the production line.

INSERT FOLDED F2 STND
7000096

Approved according to SHR_SOP 04.07.08_Att A The Signature will not be printed!	Category 2
	schustn2 - Apr 05, 2022
	gilberth - Apr 12, 2022
	sovacoeo - Apr 05, 2022
	cannons3 - Apr 13, 2022

MATRIX CODES are 5.3 mm square. Right side of matrix code is 0.87" from the right edge of page. Top of the uppermost code is 5.8549" from top edge of page. Top of the lower code is 7.3549" from the top edge of page.

PREPRESS: Matrix code must match part number:

Keep-out area for Roche material number(s) is:
0.25" (w) x 0.9" (h)
The keep out area begins:
x=17.8125"
y=6.9"
The numbers are centered from the vertical center point of 7.35"

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- Davies MJ, D'Alessio DA, Fradkin J, Kerman WN, Mathieu C, Mingrone G, Rossing P, Tsapas A, Wexler DJ, Buse JB. Management of Hyperglycemia in Type 2 Diabetes, 2018. A Consensus Report by the American Diabetes Association (ADA) and the European Association for the Study of Diabetes (EASD). *Diabetes Care*. 2018 Dec; 41(12):2669-2701.
- American Diabetes Association: 2. Classification and Diagnosis of Diabetes: Standards of Medical Care in Diabetes-2019. *Diabetes Care*, 42, (Suppl. 1): S13-S28 (2019).
- American Diabetes Association website: Diagnosing Diabetes and Learning about Prediabetes. http://www.diabetes.org/diabetes-basics/diagnosis/. Accessed April 22, 2019.
- IDF Clinical Guidelines Task Force. Global guideline for Type 2 diabetes. Brussels: International Diabetes Federation, 2012.
- Definition and diagnosis of diabetes mellitus and intermediate hyperglycemia: report of a WHO/IDF consultation. WHO, Geneva 2006 (ISBN 92 4 159493 4, ISBN 978 92 4 159493 6).
- D’Orazio et al.: “Approved IFCC Recommendation on Reporting Results for Blood Glucose (Abbreviated);” *Clinical Chemistry* 51:9 1573-1576 (2005).

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Roche Diabetes Care GmbH Sandhofer Strasse 116 68305 Mannheim, Germany www.accu-check.com